



## IFUN Dental Surgical Guides Resin

[Read More](#)

**SKU:** 3DWGGIWSFNWXX

**Price:** 949 DH

**Stock:** instock

**Categories:** [Dental Resins](#), [Resins](#)

**Tags:** [#dental](#), [3d](#), [print](#), [resin](#)

### Product Description

IFUN 3165 New R&D Dental Surgical Guides Resin for LCD DLP 3D Printer 405nm UV Cure Rapid Liquid polymer 4K Resin

#3165 Designed for 405nm wavelength LCD DLP printers to 3d print dental Surgical Guides that have higher requirements for mechanical

[caption id="attachment\_6682" align="aligncenter" width="790"]

IFUN 3165 New R&D Dental Surgical Guides Resin for LCD DLP 3D Printer 405nm UV Cure Rapid Liquid polymer 4K Resin

IFUN 3165 New R&D Dental Surgical Guides Resin for LCD DLP 3D Printer 405nm UV Cure Rapid Liquid polymer 4K Resin[/caption] [caption id="attachment\_6680" align="aligncenter" width="790"]

IFUN 3165 New R&D Dental Surgical Guides Resin for LCD DLP 3D Printer 405nm UV Cure Rapid Liquid polymer 4K Resin

IFUN 3165 New R&D Dental Surgical Guides Resin for LCD DLP 3D Printer 405nm UV Cure Rapid Liquid polymer 4K Resin[/caption] [caption id="attachment\_6681" align="aligncenter" width="790"]

IFUN 3165 New R&D Dental Surgical Guides Resin for LCD DLP 3D Printer 405nm UV Cure Rapid Liquid polymer 4K Resin

IFUN 3165 New R&D Dental Surgical Guides Resin for LCD DLP 3D Printer 405nm UV Cure Rapid Liquid polymer 4K Resin[/caption]

properties, replacing traditional of dental guide materials. High hardness & high precision, offering a good balance of properties between rigidity and functionality, which easy to insert drill sleeves,directly after printing.

eSUN



## eResin Dental Castable 1000ml

[Read More](#)

**SKU:** 3DWM264T8I9CP

**Price:** 1,299 DH

**Stock:** instock

**Categories:** [Dental Resins](#), [Functional Resins](#), [Resins](#)

**Tags:** [#dental](#), [3d](#), [print](#), [resin](#)

### Product Description

Short curing time No deformation during modling Smooth surface High precision. Good strength and toughness Excellent molding and resistance to impact. Baking at 200°C without craking.

Image not found or type unknown



Image not found or type unknown



**Castable resin for Dental Recommended Settings For Printer:**

Castable resin for Dental Recommended Settings For Printer

	Exposure Time s	Bottom Exposure Time s	Layer Height mm	Bottom Layer Height mm	Lift D
ELEGOO Saturn	4	40	0,05	5	7
ELEGOO Mars 2 Pro	4	40	0,05	5	5
ANYCUBIC Photon Mono X	1,8	60	0,05	6	8
ANYCUBIC Photon Mono 4K	4,5	45	0,05	6	6
ANYCUBIC Photon Mono SE	3	35	0,05	6	6
VOXELAB Proxima 8.9	3,5	40	0,05	5	8
PHROZEN Sonic Mighty 4K	4,2	40	0,05	6	8
CREALITY HALOT ONE	5	50	0,05	10	6
CREALITY LD006	2,8	45	0,05	5	7